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

5-Br-PAPS, a water-soluble derivative of 5-Br-PADAP, is a highly sensitive colorimetric reagent for zinc detection. This reagent forms a red chelate with zinc at pH 7.5-9.5 ($\lambda_{max}=552$ nm, $\epsilon=13.3 \times 10^4$). 5-Br-PAPS is also utilized for the determination of Cu(II), Fe(II), Ni(II), Co(II), Pt(II), Ru(II), and Rh(II). 5-Br-PAPS is utilized to determine serum levels of Zn, Cu, Fe, Ni, and Co. For serum Zn levels,

2-(5-Bromo-2-pyridylazo)-5-[N-n-propyl-N-(3-sulfopropyl)amino]phenol, disodium salt, dehydrate

Metallo Assay Zinc Assay kit is intended for the quantitative determination of zinc in biological sample by 96-well reader use. Target: zinc (Zn^{2+}) Measuring range: 4 -1,000 $\mu g/DI$ Species: all. Sample: serum, plasma, urine, saliva, tissue extract, cell lysate, biological fluid, food, beverage etc. Apparatus: UV/Vis 96-well reader or spectrophotometer. Wavelength: 560 nm (550-580nm) Procedure: 10 mi

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- **Storage conditions:** Store at 2-8°C. Don't freeze.

Expiration: 1 year from the date of manufacture. After the bottles are opened, the kit should be used in 1 month

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[Contents of this kit]

R-1 Buffer (ready to use):	24 mL
R-2 Chelate color (ready to use):	0.5 mL
Zinc standard 200 µg/dL (ready to use):	1.2 mL

Operation

1. Sample preparation

◇ Serum or Plasma

Insoluble substances in serum and plasma samples should be removed by filtration or centrifugation. EDTA-plasma cannot be used

2. Assay preparation

(1) Bring all reagents to room temperature before use.

(2) Prepare enough working Reagent (WR).

3. Assay procedure.

Procedure using microplate reader.

(1 assay sample 242µL)

? Assay

(1) Add 230 µL of Working Reagent (WR) to each well.

(2) Add 12 µL of Distilled water (Blank) / STD (Standard)/ sample into each well and incubate at room temperature for 5 min.

(3) Read the absorbance at 560 nm (main) and 700nm(sub).

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